



Multi-Channel Pipettes

- 왼손 오른손 작동 조절이 가능
- 완벽한 세척 및 건조 가능
- 정확도 및 정밀도 보장

해밀턴 SoftGrip Multi-Channel 피펫은 5 μ L ~ 50 μ L, 30 μ L ~ 300 μ L의 두가지 볼륨 범위에서 8채널과 12채널 2개의 제품을 볼륨 조절이 가능한 형태로 제공을 하며 다채널 피펫 또한 빠른 볼륨 식별을 위한 컬러 코팅 플런저 버튼과 라벨이 표기되어 있습니다.

독특하고 긴 팁 이젝터 레버와 핸들을 기준으로 360° 회전이 가능한 피펫의 몸통은 왼손 또는 오른손으로 사용을 할 때 실험자에게 최대한의 편안함을 제공해줍니다. 또한 간단하고 쉽게 작동되는 팁 이젝터 레버는 장착되어 있는 모든 팁을 한 번의 작동만으로 모두 배출시켜 줍니다.

최고의 재료로 만들어진 SoftGrip 멀티 채널 피펫은 전체를 분해하거나 재보정할 필요 없이 고압살균(Autoclaving)이 가능합니다.

그리고 Hamilton의 Universal AdvanTip™ Precision Pipette Tip과 함께 사용할 경우 정확성과 정밀도가 보장됩니다.

8-Channel Pipettes

5 μ L–300 μ L

P/N	Volume	Increments (μ L)	Color	At 10% of Pipette Volume		At 100% of Pipette Volume	
				Accuracy within $\pm$ %	Precision within + %	Accuracy within $\pm$ %	Precision within + %
55019-50	5–50 μ L x 8	0.02	● Sandstone	2.00	2.00	0.80	0.35
55019-52	30–300 μ L x 8	0.2	● Brick Red	1.20	0.70	0.50	0.30



12-Channel Pipettes

5 μ L–300 μ L

P/N	Volume	Increments (μ L)	Color	At 10% of Pipette Volume		At 100% of Pipette Volume	
				Accuracy within $\pm$ %	Precision within + %	Accuracy within $\pm$ %	Precision within + %
55019-54	5–50 μ L x 12	0.02	● Sandstone	5.00	2.00	1.00	0.50
55019-56	30–300 μ L x 12	0.2	● Brick Red	2.40	1.20	1.00	0.50

Accessories

P/N	Product	Description
34195	Multi-Channel Pipette Stand	Constructed of heavy, solvent-resistant plastic; holds one multi-channel pipette
1987-01	Calibration Labels, 120/pk	Quickly record and view calibration data with labels on your pipettes
34002	Calibration Key	For simple calibration in the field



Multi-Channel Pipette Stand

